# 504Du

## Contents

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






## English

- A. Switch on power to drive
- B. User decision to calibrate
- C. Indication of head and tubing to which pump is currently calibrated
- D. Indication of default set up. Press enter if OK, press Step enter to change set up.
- E. Decision to change setup.
- F. Calibration option by pumphead and tubing or by calibration dose.
- G. Reference Calibration flow chart.
- H. After selecting calibration setting return to Main Menu.
- I. Set direction and speed. Start to run pump.

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## English

- J. Prompt to set direction of rotation, speed
- K. Set direction and speed at keypad
- L. Press Start to start calibration dose
- M. Wait Press Stop when ready
- N. After 15 seconds prompt will change to Press stop
- O. Indication of current calibration settings.
- P. Enter actual volume/weight physically measured to allow for viscosity of fluid. Use Speed increment or decrement keys to increase or decrease the displayed volume/weight
- Q. Go to the flow chart outlining Start-up.

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## Error Messages

The pump will retain data such as pumphead, tube size, calibration, rpm and direction of rotation. If the data has been corrupted, the following two screens will appear for 2 seconds each in a continuous cycle until **Stop** or **Enter** are pressed.

If the pump detects any other fault condition, the following two screens will appear in a continuous cycle until **Stop** or **Enter** are pressed.

This screen indicates that the motor has stalled or that a loss of tachometer pulses has been selected. Clear the cause of stalling and power the pump off and on. If the error persists seek qualified assistance.

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
---

**Watson-Marlow, Bioprene and Marprene** are trademarks of **Watson-Marlow Limited**.

**Tygon** is a trademark of the **Norton Company**.

**Warning, These products are not designed for use in, and should not be used for patient connected applications.**

The information contained in this document is believed to be correct but **Watson-Marlow Limited** accepts no liability for any errors it contains, and reserves the right to alter specifications without notice.

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---

## Product Use and Decontamination Certificate

In compliance with the **UK Health & Safety at Work Act** and the **Control of Substances Hazardous to Health Regulations** you, the user are required to declare the substances which have been in contact with the product(s) you are returning to Watson-Marlow or any of its subsidiaries or distributors. Failure to do so will cause delays in servicing the product. Therefore, **please complete this form** to ensure that we have the information **before** receipt of the product(s) being returned. **A FURTHER COPY MUST BE ATTACHED TO THE OUTSIDE OF THE PACKAGING CONTAINING THE PRODUCT(S).** You, the user, are responsible for cleaning and decontaminating the product(s) before returning them.

Please complete a separate Decontamination Certificate for each pump returned.

**RGA No:** .....

1. Company .....

Address .....

..... Postcode .....

Telephone ..... Fax number .....

2. Product .....

2.1 Serial number .....

2.2 Has the product been used?

YES		NO	
-----	--	----	--

If yes, please complete all the following Sections. If no, please complete Section 5 only

3. Details of substances pumped

3.1 Chemical names

(a) .....

(b) .....

(c) .....

(d) .....

3.2 Precautions to be taken in handling these substances

(a) .....

(b) .....

(c) .....

(d) .....

3.3 Action to be taken in the event of human contact

(a) .....

(b) .....

(c) .....

(d) .....

3.4 Cleaning fluid to be used if residue of chemical is found during servicing

(a) .....

(b) .....

(c) .....

(d) .....

Note: Please describe current faults .....

.....

.....

.....

4. I hereby confirm that the only substances(s) that the equipment specified has pumped or come into contact with are those named, that the information given is correct, and the carrier has been informed if the consignment is of a hazardous nature.

5. Signed .....

Name .....

Position .....

Date .....