# Contents

EN

<b>1 Original instructions</b>	<b>4</b>
<b>2 Introduction</b>	<b>4</b>
<b>3 Pump models</b>	<b>5</b>
<b>4 Hazardous environments</b>	<b>7</b>
<b>5 Operating parameters</b>	<b>8</b>
5.1 Humidity conditions	9
5.2 Fluid temperature range	9
<b>6 Performance</b>	<b>10</b>
<b>7 Flow rates</b>	<b>11</b>
<b>8 Tube life</b>	<b>11</b>
<b>9 Potential pump hazards</b>	<b>12</b>
<b>10 Installation guidelines</b>	<b>12</b>
<b>11 Installation</b>	<b>13</b>
11.1 Installation with adapter plate (DEA1609A)	13
<b>12 Maximum number of pumpheads</b>	<b>14</b>
<b>13 Extension pumphead</b>	<b>15</b>
13.1 Fitting an extension pumphead	15
13.2 Earthplate assembly kit	16
13.3 Earthing	17
<b>14 Operation</b>	<b>18</b>
14.1 Opening the pumphead	19
<b>15 Tube loading</b>	<b>20</b>
<b>16 Dimensions</b>	<b>21</b>
16.1 Pumphead dimensions	21
16.2 Extension pumphead dimensions	21
16.3 Backplate and shaft dimensions	21
<b>17 Materials of construction</b>	<b>22</b>
<b>18 Special conditions of use for safe operation—Tube materials suitable for use with this equipment</b>	<b>23</b>
<b>19 Product codes</b>	<b>25</b>
19.1 313, 314 1.6mm wall thickness	25

---

<b>20 Maintenance, servicing and cleaning requirements</b> .....	<b>26</b>
<b>21 Summary of modifications</b> .....	<b>27</b>
<b>22 Replacements</b> .....	<b>27</b>
<b>23 ATEX marking</b> .....	<b>28</b>

Watson-Marlow Limited  
Falmouth  
Cornwall  
TR11 4RU  
England

## EU declaration of conformity

- 300 Series, ATEX compliant, peristaltic pumpheads.
- Manufacturer: Watson Marlow Ltd, Bickland Water Road, Falmouth, TR11 4RU, UK
- This declaration of conformity is issued under the sole responsibility of the manufacturer.

4. The following models and versions of covered:

313DA*	033.3411.A00	3 roller, tongued drive fitting for 1.6mm wall tubing
313D2A*	033.3511.A00	3 roller, tongued drive fitting for 2.4mm wall tubing
314DA*	033.4411.A00	4 roller, tongued drive fitting for 1.6mm wall tubing
314D2A*	033.4511.A00	4 roller, tongued drive fitting for 2.4mm wall tubing
313BA*	033.3421.A00	3 roller, bareshaft fitting for 1.6mm wall tubing
313B2A*	033.3521.A00	3 roller, bareshaft fitting for 2.4mm wall tubing
314BA*	033.4421.A00	4 roller, bareshaft fitting for 1.6mm wall tubing
314B2A*	033.4521.A00	4 roller, bareshaft fitting for 2.4mm wall tubing
313XA*	033.3431.A00	3 roller, extension head for 1.6mm wall tubing
314XA*	033.4431.A00	4 roller, extension head for 1.6mm wall tubing
313XBA*	033.3441.A00	3 roller, bareshaft extension head for 1.6mm wall tubing
314XBA*	033.4441.A00	4 roller, bareshaft extension head for 1.6mm wall tubing
313X2A	033.3531.A00	3 roller, extension head for 2.4mm wall tubing
313XB2A	033.3541.A00	3 roller, bareshaft extension head for 2.4mm wall tubing
314X2A	033.4531.A00	4 roller, extension head for 2.4mm wall tubing
314XB2A	033.4541.A00	4 roller, bareshaft extension head for 2.4mm wall tubing
313B2A	033.3521.AM0	3 roller, long bareshaft fitting for 2.4mm wall tubing
313XB2A	033.3541.AM0	3 roller, long bareshaft extension head for 2.4mm wall tubing

- The object of the declaration described above is in conformity with the relevant Union harmonisation legislation: Directive 2014/34/EU (the "ATEX" directive).  
The pumpheads carry the following ATEX categorisation:

- This declaration applies to the pump when using the Watson-Marlow tubing stated within the pumphead manual and in accordance with the special operating instructions provided in the manual. The use of any other tubing material in the pump would invalidate this declaration.
- Harmonised standards used:  
EN 80079-36:2016  
EN 80079-37:2016
- Full details of the conformity assessment procedure can be found in the technical reference file, "ATEX-WM". In accordance with the requirements of Directive 2014/34/EU a copy of this file has been archived with the following notified body:  
SGS Fimko Oy, FI-00211 Helsinki, Finland.

Signed for and behalf of:

Watson Marlow Ltd  
Falmouth, 31<sup>st</sup> October 2019

Simon Nicholson, Managing Director, Watson-Marlow Limited  
Watson-Marlow Fluid Technology Group  
Telephone +44 (0) 1326 370370  
A Spirax-Sarco Engineering plc company

## 1 Original instructions

The original instructions for this manual have been written in English. Other language versions of this manual are a translation of the original instructions.

## 2 Introduction

Directive 2014/34/EU, commonly known as the ATEX directive, carries obligations to the person who places equipment on the market, in the EU territory, for use in potentially explosive environments. A number of pumps manufactured by Watson- Marlow are suitable for use in hazardous environments; ATEX compliant pumps from the 300 series are listed below, in Section 2 "Pump models".

All of Watson-Marlow's 313/4 ATEX pumps have been rated as II 2G Ex h IIB T4 Gb X. The 313 and 314 pumphead range is intended for use in gas based environments only. The pumpheads are sold individually (for use with customer specified ATEX drives) and in defined configurations supplied complete with ATEX drives.

**This document provides specific ATEX information and should be used in conjunction with the original user manual.**

### 3 Pump models

The following 313/314 pumpheads are covered by this manual:

313 and 314 ATEX Pumpheads			
313DA*	033.3411.A00	3 roller, tongued drive fitting for 1.6mm wall tubing	
313D2A*	033.3511.A00	3 roller, tongued drive fitting for 2.4mm wall tubing	
314DA*	033.4411.A00	4 roller, tongued drive fitting for 1.6mm wall tubing	
314D2A*	033.4511.A00	4 roller, tongued drive fitting for 2.4mm wall tubing	
313BA*	033.3421.A00	3 roller, bareshaft fitting for 1.6mm wall tubing	
313B2A	033.3521.AM0	3 roller, long bareshaft fitting for 2.4mm wall tubing	
313B2A*	033.3521.A00	3 roller, bareshaft fitting for 2.4mm wall tubing	
314BA*	033.4421.A00	4 roller, bareshaft fitting for 1.6mm wall tubing	
314B2A*	033.4521.A00	4 roller, bareshaft fitting for 2.4mm wall tubing	
313XA*	033.3431.A00	3 roller, extension head for 1.6mm wall tubing	
313X2A	033.3531.A00	3 roller, extension head for 2.4mm wall tubing	
314XA*	033.4431.A00	4 roller, extension head for 1.6mm wall tubing	
314X2A	033.4531.A00	4 roller, extension head for 2.4mm wall tubing	
313XBA*	033.3441.A00	3 roller, bareshaft extension head for 1.6mm wall tubing	
313XB2A	033.3541.A00	3 roller, bareshaft extension head for 2.4mm wall tubing	
313XB2A	033.3541.AM0	3 roller, long bareshaft extension head for 2.4mm wall tubing	
314XBA*	033.4441.A00	4 roller, bareshaft extension head for 1.6mm wall tubing	
314XB2A	033.4541.A00	4 roller, bareshaft extension head for 2.4mm wall tubing	

\* NB: these pumpheads are compatible with a number of Watson-Marlow cased drives. However, the Watson-Marlow cased drives are not ATEX compliant and must not be used in hazardous locations.

**Any 313 or 314 pumpheads not listed here are NOT suitable for use in hazardous environments (at date of issue).**

## 4 Hazardous environments

All of Watson-Marlow's 313/4 ATEX pumps have been rated as II 2G Ex h IIB T4 Gb X equipment under the definitions of 2014/34/EU:

"Group II, Category 2 products should be designed to be capable of remaining within their operational parameters, as stated in the instruction manual, and based on a high level of protection for their intended use, in areas in which explosive atmospheres caused by mixtures of air and gases, vapours, mists or air/dust mixtures are likely to occur."

Furthermore, "The explosion protection relating to this category must function in such a way as to provide a sufficient level of safety even in the event of equipment with operating faults or in dangerous operating conditions which normally have to be taken into account".

Watson-Marlow pumps must not be used in the underground parts of mines, and in surface installations of such mines, likely to become endangered by firedamp and/or combustible dust.

NB: where two or more items of ATEX equipment are combined, the complete assembly shall carry the same rating as the lowest ranking individual piece of equipment.

EN

## 5 Operating parameters

The following tube materials can be used with 313 and 314 ATEX pumpheads. Bore sizes range from 1.6mm to 8.0mm, with an option of 1.6mm or 2.4mm wall thickness:

<b>Tubing: working temperature range</b>	
Marprene	5 – 80C
Bioprene	5 – 80C
Pumpsil platinum-cured silicone	-20 – 80C
STA-PURE PCS	0 – 80C
STA-PURE PFL	0 – 80C

The following parameters define the boundary of the safe working envelope - these values must not be exceeded (ATEX compliance will be invalidated):

<b>313 and 314 pumpheads</b>	
Ambient temperature range	5 – 40C
Max peak pressure (0.5 – 4.8mm bore)	2 bar
Max peak pressure (6.4 – 8.0mm bore)	1 bar
Max continuous speed, 313	400 rpm
Max continuous speed, 314	300 rpm
Corrosion resistance	See "Materials of construction" on page 22
Tube life	See "Tube life" on page 11

**WARNING! Do not run dry for excessive periods. Roller and tubing temperatures can exceed normal operating range. Tube life will be reduced increasing the chance of premature failure.**

**WARNING! Do not run pumphead against a dead-end condition (closed discharge). This can lead to excessive roller and tubing temperatures and pressures in excess of the limits in the table above. Tube life will be reduced increasing the chance of premature failure.**

When two or more items of ATEX equipment are combined, the permissible operating envelope will be determined by the narrowest range after considering all values for a given parameter.

## 5.1 Humidity conditions

The pumpheads are designed for use at non-condensing humidity levels between 80% and 35%.




## 5.2 Fluid temperature range

For all tube types, the maximum allowable fluid temperature is 65°C.

**WARNING! Do not run pumpheads at speeds higher than the rated maximum continuous speed (313 max. 400 rpm, 314 max 300 rpm). This can lead to excessive roller and tubing temperatures. Tube life will be reduced, increasing the chance of premature failure**

## 6 Performance

## 7 Flow rates

### Flow rates

(H<sub>2</sub>O, 20C, 0 bar)

Bore (mm)	313 (ml/rev)	314 (ml/rev)
1.6	0.27	0.25
3.2	1.00	0.85
4.8	2.20	1.90
6.4	3.60	3.00
8.0	5.00	4.00

## 8 Tube life

A number of factors contribute to the life of the tubing:

Factors influencing tube life

Normal tube fatigue - dependant on tube size and material.

Incorrect tube loading - see "Tube loading" on page 20 for guidance.

Excess working pressure - see "Operating parameters" on page 8.

Chemical incompatibility - a tubing compatibility guide can be found on <https://www.wmfts.com/en/support/chemical-compatibility-guide/>. Immersion kits are available from Watson-Marlow for testing.

For each application it is strongly recommended that tube life should be determined by trials, prior to installation in a hazardous environment. If this is not possible, or if there is any doubt in terms of tube life then the following hazards should be recognised before installing a pump in a potentially explosive atmosphere:

- Chemical reaction between pumped fluid and pump materials - the materials of construction are listed in "Materials of construction" on page 22.
- Pumped fluid can be ignited by surface temperature of rollers - all Watson-Marlow's ATEX equipment has been rated as T4 (meaning that even under worst-case operating conditions the maximum surface temperature will not exceed 135C).

## 9 Potential pump hazards

As part of the requirements of 2014/34/EU all potential hazards, including expected malfunctions, have been identified and subjected to a risk assessment. In order to prevent these ignition sources becoming effective a number of changes have been implemented. In addition to engineering modifications, the changes include additional operating instructions in order to specify correct usage in hazardous locations.

### Recognised ignition sources

Surface temperatures of rollers and spindles
Tube burst and subsequent spilling of pumped fluid
Mechanical failure of rotor hub
Exothermic chemical reaction
Electrostatic discharge
Bearing failure

## 10 Installation guidelines

**This symbol, used on the pumphead means: Caution, refer to accompanying documents.**

**This symbol, used on the pumphead means: Do not allow fingers to contact moving parts.**

The 313/314 pumpheads achieve ATEX II 2G Ex h IIB T4 Gb X equipment compliance for electrostatic discharge from non-conductive parts without modification. This is because the 313/314 pumphead projected surface area in any direction is less than 10,000mm<sup>2</sup>. EN ISO 80079-36:2016 par 6.7.5 for Gas Group IIB equipment states that there is a danger of electrostatic charge under normal conditions of use for equipment with greater projective area than 10,000mm<sup>2</sup>.


The bare shaft models need to be installed with the bareshaft connected to earth, this is usually achieved via the use of a suitable ATEX drive, it is the customers responsibility to ensure safe operation.

The effect of this is that a MAXIMUM OF TWO 313/314 PUMPHEADS (to suit 1.6mm tubing only) can be used in a ganged, multi-channel, Group IIB equipment configuration. If more than two pumpheads are to be ganged together then additional earthing plates will need to be installed, please refer to "Maximum number of pumpheads" on page 14 and "Extension pumphead" on page 15.


If the tubing, external to the pumphead, comes into dynamic contact with another surface or is in proximity to an electrical/magnetic field, the tubing could become electrostatically charged. Tubing should be positioned to prevent this from occurring.

## 11 Installation


EN

### 11.1 Installation with adapter plate (DEA1609A)

## 12 Maximum number of pumpheads

The table states the maximum number of 1.6mm wall thickness pumpheads that can be stacked together on a single drive shaft. Only Watson-Marlow tubing should be used in the pumphead.

<b>Marpene, Bioprene, PureWeld XL, Pumpsil</b>		Tube bore (mm) / Number of heads				
Drive	Pressure (bar)	1.6	3.2	4.8	6.4	8.0
0-350 rpm	1	5	4	2	2	1
0-350 rpm	2	5	4	2	NA	NA

Do not stack pumpheads at drive speeds that exceed 350rpm


The table states the maximum number of 2.4mm wall thickness pumpheads that can be stacked together on a single drive shaft. Only Watson-Marlow tubing should be used in the pumphead.

<b>Marpene, Bioprene, PureWeld XL, Pumpsil</b>		Tube bore (mm) / Number of heads				
Drive	Pressure (bar)	1.6	3.2	4.8	6.4	8.0
0-350 rpm	1	4	3	2	2	1
0-350 rpm	2	4	2	1	NA	NA

Do not stack pumpheads at drive speeds that exceed 350rpm

## 13 Extension pumphead

### 13.1 Fitting an extension pumphead


The coupling on a new pumphead is supplied pre-greased.

#### See also

"Maintenance, servicing and cleaning requirements" on page 26

## 13.2 Earthplate assembly kit

Item	Description	Quantity
1	Small feet	2
2	Medium feet	2
3	Large feet	2
4	Drive shaft extension	1
5	Nut	4
6	Washer	2
7	Earthing plate	1
8	Mounting plate	1
9	Earth strap	1


### 13.3 Earthing

EN

To prevent accumulation of static charge you must install an earthing plate between every two pumpheads. The earthing plates should be connected to an earth point on the drive, earth connections should meet the requirements of the ATEX directive. Continuity to earth must be tested prior to use.

- Earth plate assembly kit part number: DEAZ3087A.
- Ensure the drive shaft extension is fitted with the Earth plate kit

#### See also


"Maximum number of pumpheads" on page 14

## 14 Operation

1. Tube clamp setting

2. Tube wall thickness

3. Set for tube bore size


## 14.1 Opening the pumphead

EN


There are moving parts inside the pumphead. Before opening the pumphead guard, ensure that the following safety directions are followed:

1. Ensure that any motor drive connected to the pumphead is isolated from any electrical or compressed air supply.
2. Ensure that there is no pressure in the pipeline
3. If a tube failure has occurred, ensure that any fluid in the pumphead has been allowed to drain to a suitable vessel, container or drain
4. Ensure the pumphead is isolated from the fluid supply
5. Ensure that appropriate Personal Protective Equipment (PPE) is worn


## 15 Tube loading


## 16 Dimensions


### 16.1 Pumphead dimensions

### 16.2 Extension pumphead dimensions

### 16.3 Backplate and shaft dimensions

## 17 Materials of construction

Description	Material	Finish
Body - rear	PPA IXEF 1022	
Body - front	PPA IXEF 1022	
Track	PPA IXEF 1022	
Lever	PPA IXEF 1022	
Locking lever	Nylon 66 (30% glass fill)	
Clamp - tube	PPA IXEF 1022	
Rotor	Nylon 66 (30% glass fill)	
Roller	Nylatron GS	
Adjuster - tube clamp	Nylon 6	
Knob - tube clamp adj	Nylon 6	
Clamp - tube 313/314	PPA IXEF 1022	
Plate - adaptor	PPA IXEF 1022	
Spindle	BS970.230MO7PB (EN1A)	Nickel plate
Shaft	Stainless steel 303S31F9	
Ball bearing	Steel (1% carbon, 1.5% chromium)	
Bearing	Sintered bronze	
Spring	Stainless steel	
Pin - Groverlok	Mild steel	
Screws	Stainless steel	
Nuts	Stainless steel	
Screw - Plastite 45	Mild steel	Zinc plate

The above materials have been carefully selected and have a well proven track record. However, if there are any aggressive chemicals present then it is imperative that a risk assessment is conducted. This must not be limited to just the pumped fluid but should also include any other aggressive fluids in the intended operating environment.

## 18 Special conditions of use for safe operation—Tube materials suitable for use with this equipment

EN

The 313/314 ATEX pumpheads have been rated as

II 2G EX h IIB T4 Gb X.

The "X" denotes that users need to follow special operating instructions to achieve the ATEX rating. The special instruction in this case relates to the tubing which is permissible for use in the pump. Only the Watson-Marlow tubing listed below is to be used in these products to ensure suitability for ATEX. All the Watson-Marlow tubing listed below has been electrostatically tested in accordance with EN80079-36:2016, 6.7.5(b) and has been found to be acceptable for use in IIB gas environments. Any use of other manufacturers tubes, or Watson-Marlow tube materials/sizes not listed below is a breach of these operating instructions and may result in the equipment being unsuitable for the specified environment.

The following Watson-Marlow tube materials are suitable for use in 313/314 ATEX pumphead applications:

- Marprene
- Bioprene
- Neoprene
- PureWeld XL
- STA-PURE PCS
- STA-PURE PFL.

Watson-Marlow recommend that the length of the tube used for the application is kept to the minimum required for the pump to operate and connect to the users system. If the user requires longer lengths of tubing to be used, it is their responsibility to ensure that the system is still in compliance and suitable for the ATEX zone.

<b>Permissible tube sizes for IIB rating:</b>			
Series	Bore (mm)	Wall (mm)	OD (mm)
313,314	8.0	2.4	12.8
313,314	6.4	2.4	11.2
313,314	4.8	2.4	9.6
313,314	3.2	2.4	8.0
313,314	1.6	2.4	6.4
313,314	8.0	1.6	11.2
313,314	6.4	1.6	9.6
313,314	4.8	1.6	8.0
313,314	3.2	1.6	6.4
313,314	2.4	1.6	5.6
313,314	1.6	1.6	4.8

## 19 Product codes

### 19.1 313, 314 1.6mm wall thickness

mm	"	#	Marprene	Bioprene	Pumpsil	Pureweld XL
1.6	1/16	14	902.0016.016	933.0016.016	913.A016.016	941.0016.016
3.2	1/8	16	902.0032.016	933.0032.016	913.A032.016	941.0032.016
4.8	3/16	25	902.0048.016	933.0048.016	913.A048.016	941.0048.016
6.4	1/4	17	902.0064.016	933.0064.016	913.A064.016	941.0064.016
8.0	5/16	18	902.0080.016	933.0080.016	913.A080.016	941.0080.016

mm	"	#	GORE™ STA-PURE Series PCS	GORE™ STA-PURE Series PFL	Neoprene
1.6	1/16	14	961.0016.016	966.0016.016	920.0016.016
3.2	1/8	16	961.0032.016	966.0032.016	920.0032.016
4.8	3/16	25	961.0048.016	966.0048.016	920.0048.016
6.4	1/4	17	961.0064.016	966.0064.016	920.0064.016
8.0	5/16	18	961.0080.016	966.0080.016	920.0080.016

### 313,314 2.4mm wall thickness


mm	"	#	Marprene	Bioprene	Pumpsil	GORE™ STA- PURE Series PFL
1.6	1/16	14	902.0016.024	933.0016.024	913.A016.024	966.0016.024
3.2	1/8	16	902.0032.024	933.0032.024	913.A032.024	966.0032.024
4.8	3/16	25	902.0048.024	933.0048.024	913.A048.024	966.0048.024
6.4	1/4	17	902.0064.024	933.0064.024	913.A064.024	966.0064.024
8.0	5/16	18	902.0080.024	933.0080.024	913.A080.024	966.0080.024

## 20 Maintenance, servicing and cleaning requirements

If aggressive liquids are spilled onto the pumphead, remove any tubing from the pumphead, remove the pumphead and clean with a mild detergent.

To limit the build-up of dust (which can become electrostatically charged and/or heated by friction), remove any tubing from the pumphead, remove the pumphead and clean with a mild detergent.

Occasionally check the moving parts of the rotor for freedom of movement and lubricate pivot points and roller bearing surfaces. If the pumphead is regularly wet or washed, lightly lubricate the roller bearing surfaces. Use Teflon lubricant or FS DuraTech perfluoroether grease or similar. (Available from WMFTS in 10ml syringes part number OG0035.)

\*Grease the lever with FS DuraTech perfluoroether grease, as required after washdowns

All pumpheads should be inspected weekly for any damage, and to ensure no debris is present inside the pumphead. The rotor should be examined monthly to determine any excessive wear of the drive shaft coupling. If the rotor has backlash greater than 20 degrees then the pumphead should be replaced. The grease on the drive shaft coupling should be cleaned off and replaced every three months. WMFTS recommend Castrol Moly grease, protector series.

Ensure that regular cleaning is carried out using a damp cloth to avoid the potential for a build-up of static charge.

## 21 Summary of modifications

- ATEX features of 313/314 pumpheads
- ATEX label: This is a requirement of the Directive and includes the ATEX rating for the pumphead (II 2G IIB h T4 Gb X).
- ATEX specific manual

## 22 Replacements

Spares and replacements should be ordered through Watson-Marlow pumps or through an official representative. Only Watson-Marlow spares and replacements should be used in order to guarantee continued compliance with the ATEX directive.

Watson-Marlow's policy is to provide spare parts for all products for a minimum of 7 years from discontinuation. The ability to implement this policy is not entirely within Watson-Marlow's control and cannot be guaranteed, but every effort will be made to honour this policy.

## 23 ATEX marking

313 and 314 ATEX (A) pumpheads have been marked with these:

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Watson-Marlow Limited  
Falmouth  
Cornwall  
TR11 4RU  
England

## EU declaration of conformity

- 300 Series, ATEX compliant, peristaltic pumpheads.
- Manufacturer: Watson Marlow Ltd, Bickland Water Road, Falmouth, TR11 4RU, UK
- This declaration of conformity is issued under the sole responsibility of the manufacturer.
- The following models and versions of covered:

313DA*	033.3411.A00	3 roller, tongued drive fitting for 1.6mm wall tubing
313D2A*	033.3511.A00	3 roller, tongued drive fitting for 2.4mm wall tubing
314DA*	033.4411.A00	4 roller, tongued drive fitting for 1.6mm wall tubing
314D2A*	033.4511.A00	4 roller, tongued drive fitting for 2.4mm wall tubing
313BA*	033.3421.A00	3 roller, bareshaft fitting for 1.6mm wall tubing
313B2A*	033.3521.A00	3 roller, bareshaft fitting for 2.4mm wall tubing
314BA*	033.4421.A00	4 roller, bareshaft fitting for 1.6mm wall tubing
314B2A*	033.4521.A00	4 roller, bareshaft fitting for 2.4mm wall tubing
313XA*	033.3431.A00	3 roller, extension head for 1.6mm wall tubing
314XA*	033.4431.A00	4 roller, extension head for 1.6mm wall tubing
313XBA*	033.3441.A00	3 roller, bareshaft extension head for 1.6mm wall tubing
314XBA*	033.4441.A00	4 roller, bareshaft extension head for 1.6mm wall tubing

- The object of the declaration described above is in conformity with the relevant Union harmonisation legislation:

Directive 2014/34/EU (the "ATEX" directive).

The pumpheads carry the following ATEX categorisation:

- This declaration applies to the pump when using the Watson-Marlow tubing stated within the pumphead manual and in accordance with the special operating instructions provided in the manual. The use of any other tubing material in the pump would invalidate this declaration.
- Harmonised standards used:  
EN 80079-36:2016  
EN 80079-37:2016
- Full details of the conformity assessment procedure can be found in the technical reference file, "ATEX-WM". In accordance with the requirements of Directive 2014/34/EU a copy of this file has been archived with the following notified body:

SGS Fimko Oy, FI-00211 Helsinki, Finland.

Signed for and behalf of:  
Watson Marlow Ltd  
Falmouth, 31<sup>st</sup> October 2019

Simon Nicholson, Managing Director, Watson-Marlow Limited  
Watson-Marlow Fluid Technology Group  
Telephone +44 (0) 1326 370370  
A Spirax-Sarco Engineering plc company

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313DA*	033.3411.A00	3 roller, tongued drive fitting for 1.6mm wall tubing
313D2A*	033.3511.A00	3 roller, tongued drive fitting for 2.4mm wall tubing
314DA*	033.4411.A00	4 roller, tongued drive fitting for 1.6mm wall tubing
314D2A*	033.4511.A00	4 roller, tongued drive fitting for 2.4mm wall tubing
313BA*	033.3421.A00	3 roller, bareshaft fitting for 1.6mm wall tubing
313B2A*	033.3521.A00	3 roller, bareshaft fitting for 2.4mm wall tubing
314BA*	033.4421.A00	4 roller, bareshaft fitting for 1.6mm wall tubing
314B2A*	033.4521.A00	4 roller, bareshaft fitting for 2.4mm wall tubing
313XA*	033.3431.A00	3 roller, extension head for 1.6mm wall tubing
314XA*	033.4431.A00	4 roller, extension head for 1.6mm wall tubing
313XBA*	033.3441.A00	3 roller, bareshaft extension head for 1.6mm wall tubing
314XBA*	033.4441.A00	4 roller, bareshaft extension head for 1.6mm wall tubing
- The object of the declaration described above is in conformity with the relevant Union harmonisation legislation:

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