Watson-Marlow 505LA, 505LXA and 505LGA ATEX pumphead user manual

1 Declaration of conformity	3
2 Introduction	4
3 Pumphead types	5
3.1 505L ATEX pumphead description	5
4 Warranty	6
5 Information for returning pumps	7
6 Safety notes	7
7 Potential pump hazards	9
8 Pumphead assembly and installation	10
8.1 Installation	10
8.2 Assembly	13
8.3 Pump installation	16
9 Pumphead specification	17
10 Operating parameters	18
11 Tube life	19
12 Special conditions of use for safe operation - Tube materials suitable for use with this equipment	20
12.1 Permissible tube sizes for IIB rating:	21
12.2 Permissible double 'Y' tubing elements for IIB rating:	21
12.3 Permissible tubing for IIB rating	22
13 General operation	23
14 505L Tube loading	24
14.1 505L: Fitting double-Y tubing elements	24
14.2 505L: Fitting high flow element	25

15 Maintenance	28
16 Materials of construction	32
17 Summary of modifications	33
18 ATEX marking	33
18.1 Key	33
19 Dimensions	34
20 Replacements	34
21 Performance data	35
21.1 Pumping conditions	35
22 Patient-connected use—warning	35
23 Disclaimers	36

Original instructions

The original instructions for this manual have been written in English. Other language versions of this manual are a translation of the original instructions

1 Declaration of conformity

Watson-Marlow Limited
Falmouth
Cornwall
TR11 4RU
England

EU declaration of conformity

1. 505 Series, ATEX compliant, configured peristaltic pumpheads.
2. Manufacturer:
Watson Marlow Ltd
Bickland Water Road
Falmouth
TR11 4RU
UK
3. This declaration of conformity is issued under the sole responsibility of the manufacturer.
4. The following models and versions of covered:

505LA	ATEX Pumphead Offset Tracks
505LXA	ATEX Extension Pumphead Offset Tracks
505LGA	ATEX Gore Pumphead Offset Tracks
505LA	ATEX Pumphead Inline Tracks
505LXA	ATEX Extension Pumphead Inline Tracks
505LGA	ATEX Gore Pumphead Inline Tracks
5. The object of the declaration described above is in conformity with the relevant Union harmonisation legislation:
ATEX Directive 2014/34/EU

ATEX rating:

II 2G Ex h IIB T4 Gb X

This declaration applies to the pump when using the Watson-Marlow tubing stated within the pumphead manual and in accordance with the operating instructions provided in the manual. The use of any other tubing material in the pump would invalidate this declaration.

6. Harmonised standards used:
BS EN ISO 80079-36:2016
BS EN ISO 80079-37:2016
7. Notified body, SGS Fimko Oy (CE0598), P.O Box 30 FI-00211, Helsinki, Finland, holds a copy of the technical reference file "ATEX-WM", containing full details of the conformity assessment procedure.

Signed for and behalf of:
Watson Marlow Ltd
Falmouth, 31 October 2019

Simon Nicholson, Managing Director, Watson-Marlow Limited
Watson-Marlow Fluid Technology Group Telephone +44 (0) 1326 370370
A Spirax-Sarco Engineering plc company

2 Introduction

Directive 2014/34/EU, commonly known as the ATEX directive, carries obligations to the person who places equipment on the market, in the EU territory, for use in potentially explosive environments.

All of Watson-Marlow's ATEX pumps have been rated as II 2G Ex h IIB T4 Gb X under the definitions of 2014/34/EU:

- Equipment group II
- Equipment category 2
- Environment G
- Mechanical protection concepts EX h
- Gas group IIB
- Temperature class T4
- Equipment protection level Gb
- Special operating restrictions X (see section 12 "Special conditions of use for safe operation - Tube materials suitable for use with this equipment" on page 20)

"Equipment category 2 comprises equipment designed to be capable of functioning in conformity with the operational parameters established by the manufacturer and of ensuring a high level of protection.

Equipment in this category is intended for use in areas in which explosive atmospheres caused by gases, vapours, mists or air/dust mixtures are likely to occur occasionally. The means of protection relating to equipment in this category ensure the requisite level of protection, even in the event of frequently occurring disturbances or equipment faults which normally have to be taken into account."

Watson-Marlow pumps must not be used in the underground parts of mines, and in surface installations of such mines, likely to become endangered by firedamp and/or combustible dust.

As stated in the Directive, where two or more items of ATEX equipment are combined, the complete assembly shall carry the same rating as the lowest ranking individual piece of equipment.

All Watson-Marlow ATEX pumps covered by this manual are intended for use in gas based environments only.

If you are unsure about the meaning of this ATEX rating see section 18 "ATEX marking" on page 33 or contact your Watson-Marlow representative for advice. Watson-Marlow representatives can advise which rating and approvals products carry, but cannot evaluate nor recommend which product may be suitable for use in an end users hazardous installation. Only the end user or their qualified representative can confirm the ATEX rating of the equipment meets the requirements of their installation.

DANGER! - EXPLOSION HAZARD

Incorrectly selected ATEX equipment can cause fire or explosion.

3 Pumphead types

The following 500 series pumpheads are covered by this manual:

505LA ATEX pumpheads, offset tracks for single channel operation

505LA	053.4001.A00	505LA ATEX pumphead
505LXA	053.4011.A00	505LXA ATEX pumphead
505LGA	053.4010.A0G	505LGA ATEX pumphead

505LA ATEX pumpheads, inline tracks for two-channel operation

505LA	153.4001.N61	505LA ATEX pumphead, inline tracks
505LXA	153.4011.FB1	505LXA ATEX Extension pumphead, inline tracks
505LGA	153.4010.FB2	505LGA ATEX Gore pumphead inline tracks

These pumpheads are compatible with a number of Watson-Marlow cased drives. However, the Watson-Marlow cased drives are not ATEX compliant and must not be used in hazardous locations.

Any 500 series pumpheads not listed above are NOT suitable for use in hazardous environments (at date of issue).

3.1 505L ATEX pumphead description

The 505L ATEX pumpheads use the peristaltic principle by employing 6 occluding rollers assembled as a rotor, which occludes a specially designed tube and provides fluid flow by positive displacement as it rotates.

The pumphead is largely of coated metal construction, with the track and rotor body being anodised aluminium. The occluding rollers are of steel construction. The rotation of the occluding rollers is allowed using roller bearings supported by steel shafts. Materials of construction are detailed in this manual.

The pumphead is designed to be direct mounted or close coupled onto a suitable motor gearbox using the instructions detailed in this manual.

4 Warranty

Watson-Marlow Ltd ("Watson-Marlow") warrants this product to be free from defects in materials and workmanship for one year from the date of shipment, under normal use and service.

Watson-Marlow's sole responsibility and the customer's exclusive remedy for any claim arising out of the purchase of any product from Watson-Marlow is, at Watson-Marlow's option: repair, replacement or credit, where applicable.

Unless otherwise agreed in writing, the foregoing warranty is limited to the country in which the product is sold.

No employee, agent or representative of Watson-Marlow has the authority to bind Watson-Marlow to any warranty other than the foregoing unless in writing and signed by a director of Watson-Marlow. Watson-Marlow makes no warranty of the fitness of its products for a particular purpose.

In no event:

- i. shall the cost of the customer's exclusive remedy exceed the purchase price of the product;
- ii. shall Watson-Marlow be liable for any special, indirect, incidental, consequential, or exemplary damages, however arising, even if Watson-Marlow has been advised of the possibility of such damages.

Watson-Marlow shall not be liable for any loss, damage, or expense directly or indirectly related to or arising out of the use of its products, including damage or injury caused to other products, machinery, buildings, or property. Watson-Marlow shall not be liable for consequential damages, including without limitation, lost profits, loss of time, inconvenience, loss of product pumped, and loss of production.

This warranty does not obligate Watson-Marlow to bear any costs of removal, installation, transportation, or other charges which may arise in connection with a warranty claim.

Watson-Marlow shall not be responsible for shipping damage of returned items.

Conditions

- o Products must be returned by pre-arrangement to Watson-Marlow, or a Watson-Marlow approved service centre.
- o All repairs or modifications must have been made by Watson-Marlow Ltd, or a Watson-Marlow approved service centre or with the express permission in writing of Watson-Marlow, signed by a manager or director of Watson-Marlow.
- o Any remote control or system connections must be made in accordance to Watson-Marlow recommendations.
- o All PROFIBUS systems must be installed or certified by a PROFIBUS approved installation engineer.

Exceptions

- Consumable items including tubing and pumping elements are excluded.
- Pumphead rollers are excluded.
- Repairs or service necessitated by normal wear and tear or by lack of reasonable and proper maintenance are excluded.
- Products which, in the judgement of Watson-Marlow, have been abused, misused, or subject to malicious or accidental damage or neglect are excluded.
- Failure caused by electrical surge is excluded.
- Failure caused by incorrect or sub-standard system wiring is excluded.
- Damage by chemical attack is excluded.
- Ancillaries such as leak detectors are excluded.
- Failure caused by UV light or direct sunlight.
- Any attempt to disassemble a Watson-Marlow product will invalidate the product warranty.

Watson-Marlow reserves the right to amend these terms and conditions at any time.

5 Information for returning pumps

Before returning products, they must be thoroughly cleaned/decontaminated. The declaration confirming this should be completed and returned to us in advance of the item being shipped.

You are required to complete and return a decontamination declaration stating all fluids that have been in contact with the equipment being returned to us.

On receipt of the declaration, we will issue a Returns Authorisation Number. We reserve the right to quarantine or refuse any equipment that is not displaying a Returns Authorisation Number.

Please complete a separate decontamination declaration for each product and use the correct form that denotes the location you wish to return the equipment to.

A copy of the appropriate decontamination declaration can be downloaded from the Watson-Marlow website at www.wmftg.com/decon

If you have any queries then please contact your local Watson-Marlow representative for further assistance at www.wmftg.com/contact.

6 Safety notes

This safety information should be used in conjunction with the rest of this operating manual.

In the interests of safety, this pumphead should only be used by competent, suitably trained personnel after they have read and understood the manual and considered any hazard involved. If the pump is used in a manner not specified by Watson-Marlow Ltd, the protection provided by the pump may be impaired. Any person who is involved in the installation or maintenance of this equipment should be fully competent to carry out the work. In the UK this person should also be familiar with the Health and Safety at Work Act 1974.

This symbol, used on the pump and in the manual, means: Caution, refer to accompanying documents.

This symbol, used on the pump and in the manual, means: Do not allow fingers to contact moving parts.

This symbol, used on the pump and in the manual, means: Caution, hot surface.

This symbol, used on the pump and in the manual, means: Personal Protective Equipment (PPE) must be worn.

This pump must be used only for its intended purpose.

The pump must be accessible at all times to facilitate operation and maintenance. Access points must not be obstructed or blocked.

If hazardous fluids are to be pumped, safety procedures specific to the particular fluid and application must be put in place to protect against injury to persons.

Ensure the chemicals to be pumped are compatible with the pumphead, lubricant (where applicable), tubing, pipework and fittings to be used with the pump. Please refer to the chemical compatibility guide which can be found at: www.wmftg.com/chemical.

If you need to use the pump with any other chemical please contact Watson-Marlow to confirm compatibility.

Explosion hazard. Failure to comply may cause severe or even fatal injuries.

All work, e.g. transportation, storage, installation, connection, commissioning, servicing and maintenance must be performed in a non-explosive atmosphere.

Always check to ensure that an Exd motor gearbox is suitably rated for the hazardous zone area in which it is to be used, including ATEX, Ex and any other hazardous area legislation for the country in which it is being installed. Exd motors should only be installed by Exd qualified personnel.

Primary operator protection from rotating parts of the pump is provided by the pumphead guard. Note that pumphead guards differ, depending on the type of pumphead.

There are moving parts inside the pumphead. Before opening the pumphead guard, ensure that the following safety directions are followed:

1. Ensure that any motor drive connected to the pumphead is isolated from any electrical or compressed air supply.

2. Ensure that there is no pressure in the pipeline

3. If a tube failure has occurred, ensure that any fluid in the pumphead has been allowed to drain to a suitable vessel, container or drain

4. Ensure the pumphead is isolated from the fluid supply

5. Ensure that appropriate Personal Protective Equipment (PPE) is worn

7 Potential pump hazards

As part of the requirements of ATEX Directive 2014/34/EU all potential hazards, including expected malfunctions, have been identified and subjected to a risk assessment. In order to prevent these ignition sources becoming hazardous, a number of changes have been implemented. In addition to engineering modifications, the changes include comments in these instructions in order to specify correct usage in hazardous locations.

Recognised ignition sources of the pumphead

Surface temperatures of rollers and spindles

Tube burst and subsequent spilling of pumped fluid

Mechanical failure of rotor hub

Exothermic chemical reaction

Electrostatic discharge

Bearing failure

Spring failure

8 Pumphead assembly and installation

Preliminary checks

DANGER! - EXPLOSION HAZARD

Failure to inspect the pumphead unit for damage or check the ATEX labelling data may cause a fire or explosion.

Check the label on the pumphead to ensure that the pumphead type and the ATEX labelling conform to the planning of the plant or machine.

Check that all components are present. Inspect components for damage in transit. If anything is missing or damaged, contact your Watson-Marlow representative immediately.

8.1 Installation

Drive Selection

The pumphead should be connected to a motor with an ATEX rating which is equivalent or better than the pumphead rating II 2G Ex h IIB T4 Gb X.

DANGER! - EXPLOSION HAZARD

The drive and any attachments to the pumphead unit must also be ATEX compliant. Only motors with adequate ATEX Zone category according to type plate may be used.

Select a drive unit that can deliver at least the maximum torque required to run the pumphead in the application.

The intended maximum torque required for startup is 10Nm. The maximum torque required for continuous running is 2Nm.

CAUTION!

To avoid damage to pumphead, rotor shaft torque must not exceed 10Nm.

Pumphead mounting

The pumphead should be connected to the drive unit using the Close Coupled 505L Pumphead adaptor kit provided.


Close Coupled 505L Adaptor Kit

	Items	Qty.
1	505L Adaptor Assembly for Close Coupled Units	1
2	Screw M5 x 12mm, Pozidrive Pan Hd St.St.	4
3	Washer M5 St. St.	4

Note: Apply thread lock compound to the M5 fasteners and ensure all fasteners are tightened to 3Nm.

Pumphead mounting dimensions

505LA pumphead adapter plate

A 102mm (4.01")


B 124mm (4.88")

1 Thro. hole
C $\varnothing 45.15\text{mm} \pm 0.05\text{m}$
m (1.78" ± 0.002 ")

D 4 holes $\varnothing 5.5\text{mm}$
(0.22") equispaced
on a 58mm (2.28")
PCD as shown


E 45°

Drive Flange Mounting Dimensions

A	45°
B	4 holes M5 X 11mm (0.43") deep equispaced on a 58mm (2.28") PCD as shown
-	-
-	-
-	-
-	-
-	-

Drive Flange & Drive Shaft Dimensions


A	33mm (1.30")
B	10mm (0.39")
C	21mm (0.83")
D	3.4mm (0.13") [Width tolerance range max/min = 3.45mm to 3.35mm]
E	Ø10mm (0.39") k6 [ISO Transition (push) fit, dia. tolerance range max/min mm = 10.010 to 10.001]
F	Ø45mm (1.77") h6 [ISO Clearance (location) fit, tolerance range max/min mm = 45.000 to 44.984]
G	7mm (0.28")

8.2 Assembly

Fitting the Close Coupled Pump 505L Adaptor

- Ensure the drive/motor flange mounting threads are free from dirt and oil
- Fit the adaptor over the drive shaft and locating boss.
- Align the adaptor horizontally/vertically so the four location holes are aligned with the threaded gearbox holes.
- Secure the 505L Adaptor to the drive with the four retaining screws supplied, using a suitable No 2 Pozi drive screwdriver.

Fitting the 505LA Pumphead


- Ensure that a suitable braided earth cable has been fitted to an earth point on the selected drive. Check the drive manufacturers manual for a suitable earth connection point.
- Grease the drive shaft tongue.
- Position the pumphead centrally on the front face of the drive and manually turn the pumphead rotor until the slotted rotor shaft engages with the tongue of the drive shaft.
- Align the two lower pillar locating studs on the adaptor with their respective holes on the pumphead and the two retaining screw holes on the pumphead with the respective middle pillars on the adaptor.
- Secure the pumphead to the drive with the two M5 x 10mm retaining screws supplied using a suitable screwdriver.

- Secure the earthing cable braid to the pumphead.

- Position the adaptor cover over the 505L pumphead adaptor and push downwards and click into position.

The pumphead track assembly can now be fitted and earthing cable braided secured to the track.

Prevention and dissipation of electrostatic charge

All ATEX pumpheads include provision for the prevention and dissipation of electrostatic charge. In order to dissipate electrostatic charge effectively there must be sufficient electrical contact between the pumphead and the suitably earthed drive.

It is imperative that the 505LA pumpheads are earthed by connecting the earth terminal to earth (usually via a suitable point on the pump drive).

Check the effectiveness of any earth connection by measuring its electrical resistance. The resistance must not exceed 1 M Ω .

DANGER! - EXPLOSION HAZARD

Failure to ensure dissipation of static may cause fire or explosion.

The resistance from any point on the pumphead guard to the earth terminal is typically 25 Ohms.

Before commissioning:

- **Test the maximum resistance from any point on the guard to earth. The resistance must not exceed 1 MΩ.**
- **Test the maximum resistance from a metal point on the rotor assembly to the drive casework. The resistance must not exceed 1 MΩ.**

To fit an extension pumphead

The front of the 505L pumphead is fitted with a cover plate held in place with two screws. It must be removed before a second pumphead can be fitted.

- Lift the two levers at either side of the pumphead and remove the track.
- Remove the two screws securing the cover plate to the front of the fixed pumphead. A stubby crosspoint screwdriver may be convenient for this. Remove the cover plate.
- Grease the projecting shaft tongue at the rear of the extension pumphead and the slot at the front of the fixed pumphead. Position the extension pumphead so that the tongue and slot engage, and

align the body of the extension pumphead with the body of the first, fixed, pumphead. Fit the two supplied extension pumphead fixing screws at each side of the extension pumphead and tighten. Secure the cover plate to the extension pumphead using the two screws which secured it to the fixed pumphead.

- Position the tracks and secure them by closing the levers.

WARNING! Pumphead must only be operated with Track levers in closed position.

8.3 Pump installation

For a correctly engineered installation please ensure that the following guidelines are followed:

- **Do** not build a pump into a tight location without adequate airflow around the pump.
- **Do** keep delivery and suction tubes as short and direct as possible and follow the straightest route. Use bends of large radius: at least four times the tubing diameter. Ensure that connecting pipework and fittings are suitably rated to handle the predicted pipeline pressure. Avoid pipe reducers and lengths of smaller bore tubing than the pumphead section, particularly in pipelines on the suction side. Any valves in the pipeline (not usually needed) must not restrict the flow. Any valves in the flow line must be open when the pump is running.
- **Do** use suction and delivery pipes equal to or larger than the bore of the tube in the pumphead. When pumping viscous fluids use pipe runs with a bore several times larger than the pump tube.
- **Do** ensure that your system fluid supply and discharge pipework is suitable for the hazardous environment in which the pump is operating and doesn't allow for the accumulation of electrostatic charge.
- **Do** site the pump at or just below the level of the fluid to be pumped if possible. This will ensure flooded suction and maximum pumping efficiency.
- **Do** keep the pumphead track and all moving parts clean and free from contamination and debris.
- **Do** run at slow speed when pumping viscous fluids. Flooded suction will enhance pumping performance in all cases, particularly for materials of a viscous nature.
- **Do** limit peristaltic tubing to the length adjacent to the pumphead because peristaltic tubing is insulating. Electrostatic testing has been used to determine which Watson-Marlow tubing is suitable for use in hazardous environments. See section 12 "Special conditions of use for safe operation - Tube materials suitable for use with this equipment" on page 20 for more information. Earthed, conductive pipework should be used elsewhere in the system.
- **When using Marprene or Bioprene** continuous tubing, do re-tension the tube after the first 30 minutes of running.
- If unsure of an installation please contact your local Watson-Marlow representative for further assistance.
- **Tube selection:** The chemical compatibility lists published in Watson-Marlow publications are guides. If in doubt about the compatibility of a tube material and the duty fluid, request a Watson-Marlow tube sample card for immersion trials.

9 Pumphead specification

ATEX rating	II 2G Ex h IIB T4 Gb X
Operating temperature	5C to 40C (41F to 104F)
Storage temperature	-40C to 70C (-40F to 158F)
Humidity (non-condensing)	35% to 80%
dB rating	< 70dB (A) @ 1m

Note: Where specifications are listed in more than one operating manual, the lowest specification must be adhered to.

For further information please contact your Watson-Marlow representative.

10 Operating parameters

The following tube materials can be used with the 505LA pumpheads in double “Y” or continuous tubing configurations (two separate channels). Bore sizes range from 0.5mm to 9.6mm, with 2.4mm wall thickness. Only Watson-Marlow tubing should be used to guarantee continued compliance with the ATEX directive.

Tubing: working temperature range

Marprene	5C to 80C (41F to 176F)
Bioprene	5C to 80C (41F to 176F)
Pumpsil	-20C to 80C (-4F to 176F)
PureWeld XL	-20C to 80C (-4F to 176F)
GORE® STA-PURE® Pump Tubing – Series PCS	-20C to 80C (-4F to 176F)
GORE® STA-PURE® Pump Tubing – Series PFL	-20C to 80C (-4F to 176F)
Neoprene	0C to 80C (32F to 176F)

The following parameters define the boundary of the safe working envelope - these values must not be exceeded (ATEX compliance will be invalidated):

505LA ATEX pumpheads (all model variants)

Ambient temperature range	5C to 40C (41F to 104F)
Maximum allowable fluid temperature	65C (149F)
Max peak pressure (0.5 to 9.6 mm bore)	2 bar (29 psi)
Max continuous speed	350 rpm
Corrosion resistance	section 16 "Materials of construction" on page 32
Tube life	section 11 "Tube life" on the facing page

Please note that, although GORE® STA-PURE® Pump Tubing is rated at 4 bar in a standard 505LG pumphead, the ATEX version (505LGA) is limited to 2 bar.

WARNING! Do not run dry for excessive periods. Roller and tubing temperatures can exceed normal operating range.

Tube life will be reduced increasing the chance of premature failure.

WARNING! Do not run pumphead against a dead-end condition (closed discharge). This can lead to excessive roller and tubing temperatures and pressures in excess of the limits in the table above. Tube life will be reduced increasing the chance of premature failure.

WARNING! Do not run pumpheads at speed higher than the rated maximum speed. This can lead to excessive roller and tubing temperatures. Tube life will be reduced, increasing the chance of premature failure.

When two or more items of ATEX equipment are combined, the permissible operating envelope will be determined by the narrowest range after considering all values for a given parameter.

11 Tube life

A number of factors contribute to the life of the tubing:

Factors influencing tube life

Normal tube fatigue -dependent on tube size and material

Incorrect tube loading - see section 14 "505L Tube loading" on page 24.

Excess working pressure - see section 10 "Operating parameters" on the previous page

Chemical incompatibility - a table of tubing compatibility can be found on www.wmftg.com/chemical. Immersion kits are available from Watson-Marlow for testing.

For each application it is strongly recommended that tube life should be determined by trials, prior to installation in a hazardous environment. If this is not possible, or if there is any doubt in terms of tube life then the following hazards should be recognised before installing a pump in a potentially explosive atmosphere:

See section 16 "Materials of construction" on page 32 for information on materials of construction.

Chemical reaction between pumped fluid and pump materials the materials of construction are listed at www.wmftg.com/chemical.

Pumped fluid can be ignited by surface temperature of rollers - all Watson-Marlow's ATEX equipment has been rated as T4. (Meaning that even under worst-case operating conditions the maximum surface temperature will not exceed 135C (275F))

12 Special conditions of use for safe operation - Tube materials suitable for use with this equipment

The ATEX pumpheads have been rated as

11 2G EX h IIB T4 Gb X

The "X" denotes that users need to follow special operating instructions to achieve the ATEX rating. The special instruction in this case relates to the tubing which is permissible for use in the pump. Only the Watson-Marlow tubing listed below is to be used in these products to ensure suitability for ATEX. All the Watson-Marlow tubing listed below has been electrostatically tested in accordance with EN80079 36:2016, 6.7.5(b) and has been found to be acceptable for use in IIB gas environments. Any use of other manufacturers tubes, or Watson-Marlow tube materials/sizes not listed below is a breach of these operating instructions and may result in the equipment being unsuitable for the specified environment.

The following Watson-Marlow tube materials are suitable for use in 505 ATEX pumphead applications:

- PureWeld XL
- Marprene
- Bioprene
- Pumpsil
- GORE® STA-PURE® Pump Tubing – Series PCS
- GORE® STA-PURE® Pump Tubing – Series PFL
- Neoprene

Watson-Marlow recommend that the length of the tube used for the application is kept to the minimum required for the pump to operate and connect to the users system. If the user requires longer lengths of tubing to be used, it is their responsibility to ensure that the system is still in compliance and suitable for the ATEX zone.

12.1 Permissible tube sizes for IIB rating:

Series	Bore (mm)	Wall (mm)	OD (mm)
505	9.6	2.4	14.4
	8.0	2.4	12.8
	6.4	2.4	11.2
	4.8	2.4	9.6
	3.2	2.4	8.0
	1.6	2.4	6.4

12.2 Permissible double 'Y' tubing elements for IIB rating:

Product codes - double 'Y' tubing elements			
Tube bore (mm)	Marprene	GORE® STA-PURE® Pump Tubing – Series PFL	GORE® STA-PURE® Pump Tubing – Series PCS
1.6/2.4	902.E016.K24	960.E016.K24	965.E016.K24
3.2/2.4	902.E032.K24	960.E032.K24	965.E032.K24
4.8/2.4	902.E048.K24	960.E048.K24	965.E048.K24
6.4/2.4	902.E064.K24	960.E064.K24	965.E064.K24
8.0/2.4	902.E080.K24	960.E080.K24	965.E080.K24
9.6/2.4	902.E096.K24		

Product codes - double 'Y' tubing elements			
Tube bore (mm)	PureWeld XL	Pumpsil	Bioprene
1.6/2.4	942.E016.K24	913.AE16.024	933.E016.K24
3.2/2.4	942.E032.K24	913.AE32.024	933.E032.K24
4.8/2.4	942.E048.K24	913.AE48.024	933.E048.K24
6.4/2.4	942.E064.K24	913.AE64.024	933.E064.K24
8.0/2.4	942.E080.K24	913.AE80.024	933.E080.K24
9.6/2.4	942.E096.K24	913.AE96.024	933.E096.K24
	-	high flow 913.AH96.K24	-

12.3 Permissible tubing for IIB rating

15m Coil Product codes				
Tube bore (mm)	Marprene	Bioprene	Pumpsil	PureWeld XL
1.6/2.4	902.0016.024	933.0016.024	913.A016.024	-
3.2/2.4	902.0032.024	933.0032.024	913.A032.024	-
4.8/2.4	902.0048.024	933.0048.024	913.A048.024	-
6.4/2.4	902.0064.024	933.0064.024	913.A064.024	942.0064.024
8.0/2.4	902.0080.024	933.0080.024	913.A080.024	-
9.6/2.4	902.0096.024	933.0096.024	913.A096.024	942.0096.024

Watson-Marlow recommend that the length of the tube used for the application is kept to the minimum required for the pump to operate and connect to the users system. If the user requires longer lengths of tubing to be used, it is their responsibility to ensure that the system is still in compliance and suitable for the ATEX zone.

13 General operation

Opening the pumphead guard

- Isolate the pump from the mains power supply.
- Disconnect the tubing from the external pipework.
- Lift the two levers at either side of the pumphead and remove the track.
- Follow the reverse of the installation process above to remove the tubing from the pumphead.

Pre-load checks

- Before loading tubing, ensure that all rollers rotate freely, that the tube ports and location grooves are clean

Closing the pumphead guard and start-up

Ensure no debris, loose fastenings or other foreign objects are in pumphead before start-up. Failure to do so may result in an explosion hazard.

14 505L Tube loading

Tube loading instructions are given for the three possible methods of using a 505L pumphead.

For most dispensing, and any application where there must be minimal pulsing, the use of Watson-Marlow double-Y silicone tubing elements is recommended. These are pre-formed to the correct length and merge two out-of-phase flows into a single smooth output. These elements fit on the bayonet-style retaining pegs supplied with the pumphead.

In certain circumstances it may be useful to run two separate inlet tubes from your source to the pumphead. It is advisable to do this if the fluid is particularly viscous and the single inlet tube is found to restrict the flow rate, or if the pump is situated a long distance from the source reservoir and is running at high speed. A tube clamp set is available from Watson-Marlow (contact your local Watson-Marlow representative) to securely fix twin inlet tubes.

The 505L may also be used as a two-channel, six-roller, precision pumphead, with two entirely separate tubes. The tube clamp set contains two clamps for each suitable size of tubing to enable both the inlet and the delivery side to be secured. It should be noted that, when used in this way, the flow does not benefit from the smoothing created by combining two out-of-phase flows, and slight channel-to-channel variation may be found in the flow rates. Keeping suction and discharge tubing lengths the same will minimise this variation.

14.1 505L: Fitting double-Y tubing elements

If using a standard element of any size:

- Lift the two levers at either side of the pumphead and remove the track.
- Stretch the tubing element slightly and locate it over the mushroom pegs either side of the pumphead.
- Replace the track and secure it by closing the two levers.
- Ensure the braided earth lead is securely connected to both the track and pumphead body.

14.2 505L: Fitting high flow element

Note: When using Pumpsil 9.6mm elements, best results will be obtained from a High-flow element: part number 913.AH96.K24. The clamp is required to secure the element on the suction side of the pumphead.

If using a 9.6mm High-flow element:

- Lift the two levers at either side of the pumphead and remove the track.
- Remove the mushroom peg on the suction side of the pump by turning it counter-clockwise through one quarter turn and withdrawing it. Do not remove the peg on the delivery side.
- Lay the element across the rotor, with the smaller Y-connector over the end peg on the delivery side.
- Tube clamps are marked with the nominal bore for their intended tubing. Pull the element taut and firmly secure it with a 9.6mm clamp on the suction side of the pumphead, using the bayonet-style fixing previously used for the mushroom peg.
- Replace the track and secure it by closing the two levers.
- Ensure the braided earth lead is securely connected to both the track and pumphead body.

505L: Using twin inlet tubes

Improved flow rates may be obtained by fitting two inlet tubes and using a Y connector to obtain a smooth flow.

- Lift the two levers at either side of the pumphead and remove the track.
- Remove the mushroom peg on the suction side of the pump by turning it counter-clockwise through one quarter turn and withdrawing it. Do not remove the peg on the delivery side.
- Fit two lengths of tube to the top branches of a Y connector. Lay this assembly across the rotor, with the Y connector over the end peg on the delivery side.
- Tube clamps are available, marked with the nominal bore for their intended tubing. Pull the two lengths of tube taut and firmly secure them with a correctly-sized clamp on the suction side of the pumphead, using the bayonet-style fixing previously used for the mushroom peg. You may find it convenient to insert the tubes into the arches of the clamp before securing the clamp.
- Replace the track and secure it by closing the two levers.
- Ensure the braided earth lead is securely connected to both the track and pumphead body

505L: Using two independent tubes

This option may be used to allow the 505L to function as a twin-channel pumphead.

- Lift the two levers at either side of the pumphead and remove the track.
- Remove the mushroom peg on both sides of the pump by turning it counter-clockwise through one quarter turn and withdrawing it.
- Mark a length of 145mm on two lengths of tube up to a bore size of 8.0mm; or mark a length of 150mm on two lengths of 9.6mm bore tubing.
- Tube clamps are marked with the nominal bore for their intended tubing. Secure the tubes to the suction side of the pumphead with a correctly-sized clamp, aligning one of the marks on each tube with the inner face of the clamp. Secure the clamp by inserting the fixing bayonet into the socket previously used for the mushroom peg and rotating it clockwise one-quarter turn.
- Stretch the tubes slightly and secure them with the second clamp in the same way, with the second pair of marks aligned with the inside of the second clamp. You may find it convenient to insert the tubes into the arches of the second clamp before securing the clamp.
- Replace the track and secure it by closing the two levers.
- Ensure the braided earth lead is securely connected to both the track and pumphead body.

Note: After fitting new lengths of Marprene tubing, run the pump for about half an hour and then re-tension the tubing, making sure that the nominal lengths above are restored. If tubing lengths greater than those given are used, the tubing may move about within the pumphead and damage will occur.

15 Maintenance

Explosion hazard. Failure to comply may cause severe or even fatal injuries.

All work, e.g. transportation, storage, installation, connection, commissioning, servicing and maintenance must be performed in a non-explosive atmosphere.

Scheduled maintenance

- If aggressive liquids are spilled onto the pumphead, the head should be removed and cleaned with a mild detergent. It is the users responsibility to ensure chemical compatibility of the detergent with the pumped fluid. First, isolate the pump from the mains power supply and remove the track by lifting the levers at either side of the pumphead. Remove the tubing by stretching the element over the pegs (or remove the clamps for continuous tubing). To detach the pumphead from the drive, remove the two screws located in the backplate. The cleaning procedure above should be used to limit the build-up of dust (which can become electrostatically charged and/or heated by friction).
- Avoid exposing the bearing seals to solvents for more than one minute.
- Routine inspection of the moving parts of the pumphead rotor should be carried out to check for excessive wear. Roller bearings are sealed and do not require lubricating.
- All pumpheads should be inspected weekly for any damage, and to ensure no debris is present inside the pumphead.
- Because of the importance of dissipating any electrostatic charge the earthing lead should be regularly checked for signs of corrosion.

WARNING! Pumphead must only be operated with Track levers in closed position.

505L track adjustment


The track is set for use with 2.4mm wall tubing with bore sizes from 1.6mm to 9.6mm. Do not attempt to change the gap between the track and the top of the pumphead.

Track Assembly replacement

- Unfasten earth-bonding.

- Lift up the levers.

- Remove Track.

- The track assembly has now been removed. To replace the track assembly the reverse procedure should be applied. Always ensure the braided earth bond is connected securely to both the track and pumphead before operation.

Pumphead Removal

- Pull-off the mounting plate cover .


- Undo the lower earth-bonding screw.


- Remove Head retaining screws.

- Detach the Head.

- Remove the mounting plate.

- Fit the new Head and/or Track in reverse order.

Check the effectiveness of any earth connection by measuring its electrical resistance. The resistance must not exceed 1 M Ω .

DANGER! - EXPLOSION HAZARD

Failure to ensure dissipation of static may cause fire or explosion.

The resistance from any point on the pumphead guard to the earth terminal is typically 25 Ohms.

Before commissioning:

- **Test the maximum resistance from any point on the guard to earth. The resistance must not exceed 1 M Ω .**
- **Test the maximum resistance from a metal point on the rotor assembly to the drive casework. The resistance must not exceed 1 M Ω .**

16 Materials of construction

Description	Material	Finish
Spindle - roller (knurled)*	Stainless steel 303S31	
Rotor	Aluminium HP30 TF	Anodised
Cradle - extrusion *	Aluminium 6063	Painted
Cradle - extrusion		
Track - extrusion	Aluminium 6063	Anodised
Rod - locking	Stainless steel 303S31	
Spacer - locking rod	Nylatron GS	
Lever - locking rod	Stainless steel 303S31	
Shaft - drive	Stainless steel 303S31	
Shaft - drive (505LXA)	Stainless steel 303S31	
Plate - front	Aluminium 5083	Painted
Plate - rear	Aluminium 5083	Painted
Cover - adaptor	Aluminium 5083	Painted
Plate - end	Aluminium 5083	Painted
Tube - guide	Stainless steel 303S31	
Spring	Stainless steel	
Ball bearings - roller spindles	Carbon steel	
Ball bearings - drive shaft	Carbon steel	
Foot	Nylatron and rubber	
Screws	Stainless steel	
Washers	Nylon or stainless steel	
Side assembly (505LXA)	Aluminium 6063 and stainless 302S26	Painted
Adaptor assembly	Speedal 2011	Painted

The above materials have been carefully selected and have a well proven track record. However, if there are any aggressive chemicals present then it is imperative that a risk assessment is conducted. This must not be limited to just the pumped fluid but should also include any other aggressive fluids in the intended operating environment.

17 Summary of modifications

ATEX features of 505L pumpheads

Conductive rotor assembly	Roller end-caps have been replaced with conductive POM-C (ESD) Knurled roller pin: Cuts through anodising on rotor end-plates to reduce electrical resistance to earth
Earth strap	Earthing straps connect the track assemblies to the main body of the pumphead (and to an extension pumphead, if fitted)
ATEX label	This is a requirement of the Directive and includes the ATEX rating for the pumphead (II 2G Ex h IIB T4 Gb X)
ATEX manual	This manual has been produced specifically for this ATEX product and contains information for safe use.

18 ATEX marking


The 505LA pumpheads have been marked with the following labels:

18.1 Key

II	Equipment Group II for above ground areas (surface)
2G	Equipment Category 2G (Gas) - Zone 1
Ex h	Ignition protection labelling for mechanical devices
IIB	Group IIB – typical gas: Ethylene
T4	Temperature classification (Gas) ≤ 135°C
Gb	Group II (Gas); protection Level: High
X	Special conditions of use for safe operation - See section 12 "Special conditions of use for safe operation - Tube materials suitable for use with this equipment" on page 20

19 Dimensions

- A 185mm (7.28")
- B 124mm (4.88")
- C 51mm (2.01")
- D 105mm (4.13")
- E 48mm (1.89")
- F 34mm (1.34")
- G 130mm (5.12")
- H Allow 50mm (1.97") clearance for tube removal

20 Replacements

Spares and replacements should be ordered through Watson-Marlow pumps or through an official representative. Only Watson-Marlow spares and replacements should be used in order to guarantee continued compliance with the ATEX directive.

Watson-Marlow's policy is to provide spare parts for all products for a minimum of 7 years from discontinuation. The ability to implement this policy is not entirely within Watson-Marlow's control and cannot be guaranteed, but every effort will be made to honour this policy.

Please contact your local Watson-Marlow representative for assistance.

21 Performance data

Note: Flow rates quoted have been rounded for simplicity, but are accurate to within 5% -well within the normal tubing-tolerance variation of flow rate. They should therefore be taken as a guide. Real flow rates in any application must be determined empirically.

Note: The pumphead can be connected to a variable speed. Flowrates are linear so can be calculated for speeds less than 350rpm

505LA ml/min

Min	Max	Ratio	Speed range	1.6mm	3.2mm	4.8mm
0.1	350	3500:1	0.1 to 350 rpm	0.04-150	0.23-800	0.42-1500
Min	Max	Ratio	Speed range	6.4mm	8.0mm	9.6mm
0.1	350	3500:1	0.1 to 350 rpm	0.70-2500	0.90-3200	1.3-4375*

505LA USGPH

Min	Max	Ratio	Speed range	1.6mm	3.2mm	4.8mm
0.1	350	3500:1	0.1 to 350 rpm	0.001-2.4	0.004-12.6	0.007-23.3
Min	Max	Ratio	Speed range	6.4mm	8.0mm	9.6mm
0.1	350	3500:1	0.1 to 350 rpm	0.011-38.8	0.014-49.9	0.020-68.9*

*9.6mm bore through pumphead and to delivery; 12.7mm bore to source, using 913.AH96.K24

21.1 Pumping conditions

Pressure and viscosity

- All pressure values in this operating instruction, from which performance and life figures have been calculated relate to peak pipeline pressures.
- Although rated to 2 bar working pressure, this pump will generate in excess of 2 bar working pressure if pipeline restrictions are in place. In instances where it is critical that a working pressure of 2 bar is not exceeded, pressure relief valves should be installed in the pipeline.
- For the maximum pressure rating for this pumphead see section 10 "Operating parameters" on page 18. Users must ensure in the system design that these pressure limits cannot be exceeded.
- Impulse loses and pulsation can be minimised by ensuring one metre of smooth bore linear tubing is connected to the discharge port of the pumphead. This is especially important with viscous fluids and rigid pipework. It is the users responsibility to ensure that the pipework is suitable for the ATEX zone and application and that correct provision for dissipation of any electrostatic charge has been included.

22 Patient-connected use—warning

Warning, These products are not designed for use in, and should not be used for patient connected applications.

23 Disclaimers

The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the users responsibility to ensure product suitability for use within their application. Watson-Marlow, LoadSure, Qdos, ReNu, LaserTraceability, Pumpsil, PureWeld XL, Bioprene, Marprene are registered trademarks of Watson-Marlow Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB.

GORE and STA-PURE are registered trademarks of W.L. Gore and Associates.

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Watson-Marlow Limited
Falmouth
Cornwall
TR11 4RU
England

EU declaration of conformity

- 505 Series, ATEX compliant, configured peristaltic pumpheads.
- Manufacturer:
Watson Marlow Ltd
Bickland Water Road
Falmouth
TR11 4RU
UK
- This declaration of conformity is issued under the sole responsibility of the manufacturer.
- The following models and versions of covered:

505LA	ATEX Pumphead Offset Tracks
505LXA	ATEX Extension Pumphead Offset Tracks
505LGA	ATEX Gore Pumphead Offset Tracks
505LA	ATEX Pumphead Inline Tracks
505LXA	ATEX Extension Pumphead Inline Tracks
505LGA	ATEX Gore Pumphead Inline Tracks
- The object of the declaration described above is in conformity with the relevant Union harmonisation legislation:
ATEX Directive 2014/34/EU

ATEX rating:

II 2G Ex h IIB T4 Gb X

This declaration applies to the pump when using the Watson-Marlow tubing stated within the pumphead manual and in accordance with the operating instructions provided in the manual. The use of any other tubing material in the pump would invalidate this declaration.

- Harmonised standards used:
BS EN ISO 80079-36:2016
BS EN ISO 80079-37:2016
- Notified body, SGS Fimko Oy (CE0598), P.O Box 30 FI-00211, Helsinki, Finland, holds a copy of the technical reference file "ATEX-WM", containing full details of the conformity assessment procedure.

Signed for and behalf of:
Watson Marlow Ltd
Falmouth, 31 October 2019

Simon Nicholson, Managing Director, Watson-Marlow Limited
Watson-Marlow Fluid Technology Group Telephone +44 (0) 1326 370370
A Spirax-Sarco Engineering plc company

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