



Improving your PERFORMANCE

401U/DM2 e 401U/DM3 Pompe ultracompatte a controllo manuale/automatico

Caratteristiche


- 401U/DM2 and 401U/DM3 provide two and three channel flows up to 36ml/min per channel
- Highly-compact with single control potentiometer for direction, speed and start/stop
- Analogue speed control and remote switching via a 15D connector
- Reversible for easy fluid recovery
- Uses standard double segment tubing
- Adjustable occlusion (pressure)



Portate per canale (ml/min)

Codice colore	Arancio/ Nero	Arancio/ Rosso	Arancio/ Blu	Arancio/ Verde	Arancio/ Giallo	Arancio/ Bianco	Nero/ Nero	Arancio/ Arancio	Bianco/ Bianco	Rosso/ Rosso	Grigio/ Grigio	Giallo/ Giallo	Giallo/ Blu
Daimetro mm	0.13	0.19	0.25	0.38	0.50	0.63	0.76	0.88	1.02	1.14	1.29	1.42	1.52
1.2 to 12rpm	0.001-0.01	0.002-0.02	0.004-0.04	0.01-0.1	0.02-0.2	0.03-0.3	0.04-0.4	0.05-0.5	0.07-0.7	0.08-0.8	0.11-1.1	0.13-1.3	0.15-1.5
4.0 to 40rpm	0.004-0.04	0.008-0.08	0.014-0.14	0.03-0.3	0.06-0.6	0.09-0.9	0.12-1.2	0.17-1.7	0.22-2.2	0.28-2.8	0.35-3.5	0.43-4.3	0.49-4.9
10 to 100rpm	0.009-0.09	0.020-0.20	0.035-0.35	0.08-0.8	0.14-1.4	0.22-2.2	0.31-3.1	0.43-4.3	0.55-5.5	0.70-7.0	0.89-8.9	1.1-11	1.2-12

Tube code	Blu/ Blu	Verde/ Verde	Viola/ Viola	Viola/ Nero	Viola/ Arancio	Viola/ Bianco
Daimetro mm	1.65	1.85	2.05	2.29	2.54	2.79
1.2 to 12rpm	0.17-1.7	0.21-2.1	0.26-2.6	0.31-3.1	0.37-3.7	0.44-4.4
4.0 to 40rpm	0.56-5.6	0.70-7.0	0.85-8.5	1.0-10	1.2-12	1.5-15
10 to 100rpm	1.4-14	1.8-18	2.1-21	2.6-26	3.1-31	3.6-36



401U/DM2 e 401U/DM3 Pompe ultracompatte a controllo manuale/automatico

Codici	rpm	Descrizione
040.181D.M20	1.2-12	401U/DM2 (15VAC) 220-240V 50/60Hz 1ph one or two channel manifold tubing pump
040.1H1D.M20	4.0-40	401U/DM2 (15VAC) 220-240V 50/60Hz 1ph one or two channel manifold tubing pump
040.1P1D.M20	10-100	401U/DM2 (15VAC) 220-240V 50/60Hz 1ph one or two channel manifold tubing pump
040.181D.M30	1.2-12	401U/DM3 (15VAC) 220-240V 50/60Hz 1ph three channel manifold tubing pump
040.1H1D.M30	4.0-40	401U/DM3 (15VAC) 220-240V 50/60Hz 1ph three channel manifold tubing pump
040.1P1D.M30	10-100	401U/DM3 (15VAC) 220-240V 50/60Hz 1ph three channel manifold tubing pump
049.0011.000 049.0031.000 049.0051.000 049.0071.000	400AC1 400AC2 400AC3 400AC4	0-10V CW/STOP/CCW 1.5m cavo controllo remoto 4-20mA CW/STOP/CCW 1.5m cavo controllo remoto CW/STOP/CCW 1.5m cavo controllo remoto START/STOP 3m cavo controllo remoto con pedale incluso

Tubo				
Tubo Manifold per teste 410U/DM2 and 401U/DM3				
Codice colore	Dia. (mm)	Marpren* [*]	PVC	Silicone
Arancio/Nero	0.13		981.0013.000	
Arancio/Rosso	0.19		981.0019.000	
Arancio/Blu	0.25	979.0025.000	981.0025.000	
Arancio/Verde	0.38	979.0038.000	981.0038.000	
Arancio/Giallo	0.50	979.0050.000	981.0050.000	
Arancio/Bianco	0.63	979.0063.000	981.0063.000	983.0063.000
Nero/Nero	0.76	979.0076.000	981.0076.000	983.0076.000
Arancio/Arancio	0.88	979.0088.000	981.0088.000	983.0088.000
Bianco/Bianco	1.02	979.0102.000	981.0102.000	983.0102.000
Rosso/Rosso	1.14	979.0114.000	981.0114.000	983.0114.000
Grigio/Grigio	1.29	979.0129.000	981.0129.000	983.0129.000
Giallo/Giallo	1.42	979.0142.000	981.0142.000	983.0142.000
Giallo/Blu	1.52	979.0152.000	981.0152.000	983.0152.000
Blu/Blu	1.65	979.0165.000	981.0165.000	983.0165.000
Verde/Verde	1.85	979.0185.000	981.0185.000	983.0185.000
Viola/Viola	2.05	979.0205.000	981.0205.000	983.0205.000
Viola/Nero	2.29	979.0229.000	981.0229.000	983.0229.000
Viola/Arancio	2.54	979.0254.000	981.0254.000	983.0254.000
Viola/Bianco	2.79	979.0279.000	981.0279.000	983.0279.000

* For autoclavable Marprene tubing, please replace last "0" with "*" - for example 979.0238.00*
Please state the product code when ordering pumps and tubing

Specifiche

Peso (escl. Trasf.)	1.0kg
Temperatura lavoro	5°C - 40°C
Massima rumorosità	<70dBA a 1m
Control ratio	10:1
Normative	CE, BS0800, IEC335-1, EN60529 (IP21)
Alimentazione	15VAC tramite trasformatore con adattatore AC/DC

Testa e Drive

Drive	Extruded anodised aluminium
Front support plate	Anodised aluminium
Pumphead track	Extruded anodised aluminium
Pumphead rollers	Acetal
Pumphead rotor	Aluminium



Watson-Marlow-Bredel Pumps
Falmouth
Cornwall
England
TR11 4RU
Tel: +44 (0) 1326 370370
Fax: +44 (0) 1326 376009
www.watson-marlow.com
Email support4gb@watson-marlow.co.uk