

**Y-connectors - 10 pcs - PVDF/Kynar -
with material certificate
- validation pack available**


Flexicon
Liquid Filling

Sales manual

Rev.: 07

11.04.2011

AEF



Part number	A (OD) (mm)	B (ID) (mm)	For ID tube (mm)
84-012-002	2.4	1.2	0.5 – 0.8
84-012-004	3.6	1.8	1.2 – 1.6
84-012-006	5.9	3.0	3.0 – 3.2
84-012-008	7.2	3.6	3.2 – 4.8
84-012-010	9.5	4.8	4.8 – 6.5

- Material: PVDF/Kynar
- Compliance: FDA USP Class IV, FDA CFR 177.2510, BSE/TSE Free
- Sterilisation possible by: Gamma irradiation, Ethylene Oxide, Autoclave
- Validation pack: Available upon request and by signing confidentiality agreement