

**Tube for
PF6 / PD12 / FP50 / FPC50**




Sales manual

Rev.: 07

04.04.2013

CHA



Part number	Description	A (ID) (mm)	B (wall thickness) (mm)	C (length) (m)
84-103-005	ACCUSIL™ Platinum cured silicone tubing	0.5	1.6	10
84-103-008	ACCUSIL™ Platinum cured silicone tubing	0.8	1.6	10
84-103-012	ACCUSIL™ Platinum cured silicone tubing	1.2	1.6	10
84-103-016	ACCUSIL™ Platinum cured silicone tubing	1.6	1.6	10
84-103-032	ACCUSIL™ Platinum cured silicone tubing	3.2	1.8	10
84-103-048	ACCUSIL™ Platinum cured silicone tubing	4.8	2.0	10
84-103-060	ACCUSIL™ Platinum cured silicone tubing	6.0	2.1	10
84-103-080	ACCUSIL™ Platinum cured silicone tubing	8.0	2.2	10
84-104-005	ACCUSIL™ on a reel – 150 m	0.5	1.6	150
84-104-008	ACCUSIL™ on a reel – 150 m	0.8	1.6	150
84-104-012	ACCUSIL™ on a reel – 150 m	1.2	1.6	150
84-104-016	ACCUSIL™ on a reel – 150 m	1.6	1.6	150
84-104-032	ACCUSIL™ on a reel – 150 m	3.2	1.8	150
84-104-048	ACCUSIL™ on a reel – 125 m	4.8	2.0	125
84-104-060	ACCUSIL™ on a reel – 90 m	6.0	2.1	90
84-104-080	ACCUSIL™ on a reel – 65 m	8.0	2.2	65
84-311-008	Silicone tubing	0.8	1.6	25
84-311-016	Silicone tubing	1.6	1.6	25
84-311-032	Silicone tubing	3.2	1.6	25
84-311-048	Silicone tubing	4.8	2.0	25
84-311-060	Silicone tubing	6.0	2.0	25
84-311-080	Silicone tubing	8.0	2.0	25

- ACCUSIL is recommended for optimum filling accuracy
- A validation guide can be provided free of charge
- The tubing is supplied double bagged and with a material certificate
- A short delivery time is guaranteed