

COMPLIANCE SUMMARY GUIDE

Accusil™ Issue 1



1. Regulation compliance statements

1a. Gluten

Gluten is not used in the raw material or manufacture of Accusil. As such it would not be expected to be present in Accusil tubing or extracted solutions from Accusil.

1b. Animal Derivative Content Free (ADCF) – Tallow/BSE/TSE

Watson Marlow works closely with our raw materials suppliers to ensure traceability of our raw materials and production process. We can confirm that Accusil and its raw materials are manufactured from ADCF components. As such Accusil is considered Animal Derived Component Free.

1c. Melamine

Melamine is not used in the raw material or manufacture of Accusil. As such it would not be expected to be present in Accusil tubing or extracted solutions from Accusil.

1d. Phthalates

Phthalates are not used in the raw material or manufacture of Accusil. As such it would not be expected to be present in Accusil tubing or extracted solutions from Accusil.

1e. Bisphenol A (BPA)

BPA is not use in the raw material or manufacture of Accusil. As such it would not be expected to be present in Accusil tubing or extracted solutions from Accusil.

1f. Latex

Latex is not used in the raw material or manufacture of Accusil. This includes contact materials during the manufacturing process.

1g. Presence of Heavy metals

Lead, Mercury, Cadmium and Arsenic are not used in the raw material or manufacture of Accusil.

Extraction solutions obtained after prolonged contact (7 days at 40C) with Accusil were analysed using Inductively Couple Plasma – Mass Spectrometry (ICP-MS) and have shown these elements were not present.

1h. Restriction of Hazardous Substances (RoHS) EU Directives 2002/95/EC

No hazardous substances as defined by the RoHS directive have been added in the manufacture or the raw materials of Accusil.

1i. Registration, Evaluation, Authorisation of Chemicals Hazardous to Health (REACH)

Watson-Marlow's products comply with the REACH regulations. Furthermore no substances are intentionally released as a result of using our products for their intended purpose.

Accusil does not contain any Substance of Very High Concern (SVHC) as established by the European Commission HA Candidate list posted in December 2012.