



Y-connector with certificate – Part number	A (OD) (mm)	B (ID) (mm)	For ID tube (mm)
84-011-002	2.4	1.2	0.5 – 0.8
84-011-004	3.6	1.8	1.2 – 1.6
84-011-006	5.9	3.0	3.0 – 3.2
84-011-008	7.2	3.6	3.2 – 4.8
84-011-010	9.5	4.8	4.8 – 6.5
84-011-012	12.4	7.5	8.0 – 10.0
84-011-014	15.9	10.2	10.0 – 12.5

- When choosing the Y-connector, the ID should be similar to the ID of the tubing
- Material: Polypropylene
- The above product numbers include 10 Y-connectors packed in bags and a 2.1 certificate