

PF22/QC22

PF22 peristaltic filling machine


Features and benefits

- Fill volumes from 5-5000 ml with $\pm 0.5\%$ accuracy
- Changing the fluid path can be achieved in less than one minute
- Consistent filling accuracy prevents costly overfilling
- No foaming, splashing or dripping between fills
- No cross-contamination



PF22/QC22 performance

PF22 – capacity diagram (based on water)




Technical specifications

	PF22/QC22
Filling volume range	5-5000 ml
Accuracy	±0.5 %
Ingress Protection	IP31
Power supply	110/230 V AC, 50/60 Hz
Power consumption	Max 350 W
Weight	25 kg
Weight	55 lbs
Accessories	Filling nozzles, tubes, filling stand, pneumatic diving nozzle system, foot switch, printer
Compatible tubing size	3.0 to 12.5 mm ID

Materials of construction

	PF22/QC22
Casework	Anodised aluminium, Stainless steel AISI 304
Pumphead body assembly	Anodised aluminium
Pumphead roller assembly	Anodised aluminium, Stainless steel 304

PF22/QC22 dimensions



Dimensions in mm/inch

Disclaimer: The information in this document is believed to be correct, but Watson-Marlow Flexicon A/S accepts no liability for any errors it contains, and reserves the right to alter specifications without notice. WARNING: These products are not designed for use in, and should not be used for, patient-connected applications. WMArchitect and Accusil are registered trademarks.

wmfts.com/global



05 July 2024